# **ELATION™**

BALLOON DILATOR

Wireguided Balloon Dilation Catheter

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## INSTRUCTIONS FOR USE

### CONTENTS:

- Wireguided Balloon Dilation Catheter
- 0.035" (0.89mm) Guidewire
- Guidewire Locking Device

### DESCRIPTION:

The Elation™ Wireguided Balloon Dilation Catheter is a multistage balloon that provides three (3) distinct diameters at three corresponding pressures. These pressures are clearly indicated on the device packaging and are also found on an information tag attached to the balloon catheter body. Inflation of the balloon is performed by attaching an appropriate inflation system with pressure monitoring gauge (such as the BIG60 Inflation Device), to the balloon luer (Conical 6% Female Luer Lock) on the proximal portion of the balloon catheter. The Elation Wireguided Balloon Dilation Catheter can also be passed through a minimum 2.8mm working channel endoscope. Two fluoroscopic markers are located at the proximal and distal shoulders of the balloon to aid in placement of the balloon in relation to anatomical landmarks. 1 cm black striped reference marks are located on the proximal end of the catheter with a double black stripe located every 5 cm. The distal catheter has a solid 5 cm long mark to indicate that the balloon has fully exited the endoscope. This catheter comes packaged with a 0.035" (0.89 mm), floppy tip guidewire, preloaded in the guidewire lumen. The guidewire is 33 cm longer than the balloon catheter, with the excess length extending from the guidewire luer of the catheter.

### 180cm CATHETER

Ref #	Balloon OD Ø		Inflation Pressure	
	mm	F	ATM	kPa
EPB6	6-7-8	18-21-24	6-8.5-10	608-861-1013
EPB8	8-9-10	24-27-30	4-6-8	405-608-811
EPB10	10-11-12	30-33-36	4.5-6.5-8	456-659-811
EPB12	12-13.5-15	36-40.5-45	4-6.5-8	405-659-811
EPB15	15-16.5-18	45-49.5-54	4-5.5-7	405-557-709
EB18	18-19-20	54-57-60	3.5-5-6	355-507-608

### 240cm CATHETER

Ref #	Balloon OD Ø		Inflation Pressure	
	mm	F	ATM	kPa
EPCB6	6-7-8	18-21-24	6-8.5-10	608-861-1013
EPCB8	8-9-10	24-27-30	4-6-8	405-608-811
EPCB10	10-11-12	30-33-36	4.5-6.5-8	456-659-811
EPCB12	12-13.5-15	36-40.5-45	4-6.5-8	405-659-811
EPCB15	15-16.5-18	45-49.5-54	4-5.5-7	405-557-709
ECB18	18-19-20	54-57-60	3.5-5-6	355-507-608

This product contains no detectable latex.

### INDICATIONS FOR USE:

The Elation Wireguided Balloon Dilation Catheters are intended for use in adult and adolescent populations to endoscopically dilate strictures of the alimentary tract. Also indicated in adults for endoscopic dilation of the Sphincter of Oddi with or without prior sphincterotomy.

### CONTRAINDICATIONS:

None known.

### WARNING:

#### DO NOT ATTEMPT TO REPAIR.

- Never use gas or air to inflate the Elation Wireguided Balloon Dilation Catheter.
- Check for proper position of the balloon catheter using endoscopic visualization. Balloon inflation in an improper location may lead to patient injury.
- Clinicians performing fluoroscopic-guided procedures should be trained in safety measures and aware of the potential for serious radiation-induced injury caused by long periods of fluoroscopy especially in adolescent population.

### PRECAUTIONS:

- Inspect the Elation Wireguided Balloon Dilation Catheter and packaging for damage prior to use. Do not use product if opened or damaged. Confirm that the device is consistent with the package label. Contact Customer Service to report and replace damaged product.
- For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or reesterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or reesterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- The Elation Wireguided Balloon Dilation Catheter is designed to pass over a 0.035" (0.89mm) guidewire, through a minimum 2.8mm working channel of the endoscope.
- Any use for procedures other than those indicated in these instructions is not recommended.
- The safety and effectiveness of the use of this device to dilate regions of the biliary tree other than the Sphincter of Oddi have not been established.

### INSTRUCTIONS FOR USE:

1. Choose the appropriate Elation Wireguided Balloon: The recommended application is printed on the package label. When choosing diameter, the maximum diameter should not exceed the diameter of the healthy lumen. Healthy lumen diameter can be assessed endoscopically with direct visualization or via reconstructed CT scan imagery. For dilation of the Sphincter of Oddi, the balloon size should be selected based on the size of the stone and must not exceed the diameter of the distal bile duct.
2. Prepare a 60ml inflation system, with pressure monitoring gauge (such as the BIG60® Inflation Device), with sterile water, sterile saline or a contrast mixture (up to 50% contrast medium), depending on technique.

**NOTE:** Refer to the manufacturer's directions accompanying the inflation system for instructions on preparation and use.

3. Remove the Elation Wireguided Balloon Dilation Catheter from the packaging. Unclip the catheter from the connected clip located at the proximal end of the catheter. Attach the Elation Wireguided Balloon

Dilation Catheter luer to the appropriately prepped inflation system. A stopcock may be used between the connection of the catheter luer and inflation system to aid in removal of air from system.

4. To minimize balloon profile, apply vacuum to the catheter before removing the protective sheath. Remove the protective sheath and inspect the catheter for any signs of damage.

**CAUTION:** Balloon dilation catheters should be used by or under the supervision of physicians thoroughly trained in alimentary tract balloon dilation. A thorough understanding of the technical principals, clinician application and risks associated with balloon dilation of the alimentary tract is necessary before using this device.

**NOTE:** Elation Wireguided Balloon Dilation Catheter balloons have been tested to withstand 10 successive, 30 second inflations each, to maximum labeled inflation pressure.

**CAUTION:** Avoid inflating or pre-dilating the balloon catheter prior to insertion into the endoscope. If resistance is met during the procedure, do not advance the catheter without first determining the cause of the resistance and taking remedial action.

5. Maintain vacuum to the catheter during insertion through the endoscope. Advance the catheter into the endoscope using short, deliberate, 2-3 cm movements. Due to variations in endoscope construction, resistance may be experienced immediately upon entering the endoscope and again 2-3 cm before exiting the distal end of the working channel.
6. Once the balloon has exited the distal end of the endoscope and is within endoscopic view, the guidewire may be advanced beyond the distal end of the catheter.

To use the guidewire as a catheter guide:

- Confirm position of switch on the locking device is in the "ON" position.
- Advance the guidewire into the desired position beyond the catheter tip (fluoroscopy is recommended if advancing without direct visualization).
- Length of guidewire advancement can be approximated using the white marking bands on the segment of the guidewire extending from the locking device (marking bands are placed 5 cm apart).
- Advance the catheter over extended portion of the guidewire until the balloon segment is in the desired position.
- To aid in proper placement of the balloon across the stricture, match the radiopaque markers, located at the distal and proximal shoulders of the balloon with the distal and proximal ends of the stricture.

For dilation of the Sphincter of Oddi or in very tight or tortuous areas of the anatomy, the following technique may be useful:

7. If desired, the guidewire may be removed from the catheter prior to endoscope insertion. A standard 0.035" (0.89mm) guidewire may be placed through the working channel of the endoscope and across the stricture. The catheter may then be advanced over the guidewire.
8. When balloon and guidewire are in the desired location, the guidewire can be secured with the Locking Device by moving the slide to the "OFF" position.
9. Position the balloon in the appropriate location to dilate the stricture.

**PRECAUTION:** Endoscopic or Fluoroscopic visualization should be used to confirm proper placement of the catheter. Verify that the shaft segment of the catheter is within endoscopic view. This ensures that the balloon has exited the endoscope completely prior to inflation. Fluoroscopy may also be used to confirm balloon placement. Radiopaque markers are included at the proximal and distal shoulders of the balloon.

10. Balloon must be filled with fluid. After releasing vacuum from insertion, confirm air has been removed from inflation system.
11. Once the balloon is positioned across the stricture, inflate the Elation Wireguided Balloon Dilation Catheter using the prepped inflation system, to the first of the three diameter stages. The diameter of the balloon is dictated by a corresponding pressure, as indicated on the product label and information tag attached to the catheter body.
12. Watch the waist of the balloon to determine complete balloon expansion. Once inflated, maintain position for 10 to 60 seconds until desired effect is reached. Monitor pressure by utilizing the pressure gauge of the inflation system attached to the catheter balloon luer.
13. To achieve larger diameters, continue applying pressure until the remaining diameters of the Elation Wireguided Balloon Dilation Catheter have been reached. Do not over inflate past the maximum pressure indicated on the product labeling.

**WARNING:** Do not exceed the diameter of the distal bile duct when dilating the Sphincter of Oddi.

**NOTE:** Fluctuations in pressure may be observed during dilation.

These fluctuations may require additional pressure adjustments from the inflation system (a slight drop in pressure at each diameter is normal).

**WARNING:** To prevent balloon burst, do not exceed the inflation pressure given for the largest diameter on the catheter and package label. If the balloon does rupture or a significant loss of pressure within the balloon occurs, deflate the balloon completely and carefully remove the balloon and endoscope together as a unit. Do not attempt to withdraw a ruptured balloon through the endoscope. Continue the procedure with a new catheter.

#### **CATHETER REMOVAL:**

14. Using the inflation system, create a negative pressure to completely deflate the Elation Wireguided Balloon Dilation Catheter prior to removal. Confirm that the balloon has been completely deflated (approximately 10-30 seconds depending on balloon size and inflation medium) using fluoroscopic and/or endoscopic visualization.
15. Remove the Elation Wireguided Balloon Dilation Catheter.

**WARNING: THE BALLOON MUST BE THOROUGHLY DEFLATED AND ALL FLUID REMOVED PRIOR TO WITHDRAWAL.**

**PRECAUTION:** Do not pull back on the catheter until the balloon is deflated completely. For improved withdrawal, straighten the distal end of the endoscope as much as possible. Any excess bend in the working channel will increase the force needed to withdraw the Elation Wireguided Balloon Dilation Catheter through the endoscope. If excessive resistance is felt, remove the endoscope & deflated balloon catheter together as a complete unit to prevent damage to body tissue, the catheter or endoscope.

#### **DEVICE DISPOSAL:**

After use, the sheath, catheter, guidewire, locking device, inflation device & stopcock should be disposed of in a

manner consistent with standard protocols for biohazard waste disposal.

**NECESSARY ACCESSORIES:**

- Inflation system with pressure monitoring gauge – 60ml, 0-12 ATM (max. 5% error), with conical 6% male luer lock.

**INCLUDED OPTIONAL ACCESSORIES:**

- Guidewire Locking Device
- 0.035" (0.89mm) Guidewire

**STORAGE:**

Store in a cool, dry place.

**COMPLICATIONS:**

Possible complications that may result from an alimentary tract dilation procedure include, but may not be limited to:

- Perforation
- Hemorrhage
- Hematoma
- Sepsis/Infection
- Allergic reaction to contrast medium
- Pain and tenderness

**WARRANTY:**

The manufacturer warrants that reasonable care has been used in the design and manufacture of this device. This warranty is exclusive and manufacturer makes no other representations or warranties of any kind to customers, its end users, or to any third parties with respect to the device and hereby expressly disclaims any and all other warranties, express or implied, statutory or otherwise, including, but not limited to, infringement and the implied warranties of merchantability and fitness for a particular purpose, even if manufacturer is aware of such purpose. Handling and storage of this device, as well as other factors relating to the patient, diagnosis, treatment, implant procedures, and other matters beyond the control of the manufacturer, directly affect the device and the results obtained from its use. The manufacturer's obligation under this warranty is limited to the replacement of the device. Under no circumstances shall manufacturer be liable to customer or any other person or entity for any punitive, special, incidental or consequential damages directly or indirectly arising from the use of this device. The manufacturer neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. This warranty shall not apply, and manufacturer assumes no liability with respect to, devices that have been (i) modified, changed, altered, misused, mishandled, repaired, reused, reprocessed, refurbished or resterilized; (ii) subjected to improper maintenance, testing or storage, accident, tampering, or inadequate protection against shock, vibration, excessively high or low temperatures, overpressure, or physical, environmental or electrical stress; (iii) been used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the device specifications, or in an application or environment for which such device was not designed or contemplated; or (iv) distributed or used contrary to applicable federal, state, local or regulatory standards.

	Single Use
	CAUTION: Consult accompanying documents
	Maximum Inflation Pressure
	Sterilized using ethylene oxide
	RX Only: CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
	Do not use if package is damaged
	Minimum working channel

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Merit Medical Systems, Inc.
1600 West Merit Parkway, South Jordan, Utah 84095 U.S.A. 1-801-253-1600
U.S.A. Customer Service 1-800-356-3748
Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland