

TRAPCAPtain™

Specimen Collection Trap

Bronchoscope securement clip keeps the bottle from flipping to a horizontal position

Angled suction tubing allows the bottle to naturally hang in a vertical position from the bronchoscope

Suction tubing securement clip adds additional stability



MERITMEDICAL
ENDOTEK®

An easy to manage bronchial wash and BAL trap with features that prevent accidental specimen loss



TRAP CAPtain™

Specimen Collection Trap

Catalog
Number

Bottle Size

Packaged
(sterile)

TC80

80 mL

10/box

Includes specimen container and lid for transporting liquid specimens

Also available with BAL Convenience Kit



Merit Medical Endotek

1600 West Merit Parkway
South Jordan, Utah 84095
1-801-253-1600
1-800-35-MERIT
merit.com/endotek

Merit Medical Europe, Middle East, and Africa (EMEA)

Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
Tel: +31 43 358 82 22

Merit Medical Ireland Ltd.

Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Free phone for specific country

Austria	0800 295 374	Luxembourg	8002 25 22
Belgium	0800 72 906 (Dutch)	Netherlands	8000 022 81 84
	0800 73 172 (French)	Norway	800 11629
Denmark	80 88 00 24	Sweden	020 792 445
France	0800 91 60 30	Switzerland	800 9000 5060
Finland	0800 770 586	UK	8000 973 115
Germany	0800 182 0871		
Ireland (Republic)	1800 553 163		
Italy	800 897 005		



EXPANDING THE POSSIBILITIES™

401494001/A 2016-03-24