

# Segue™

Multipurpose TransProtection System™



EXPANDING THE POSSIBILITIES™



**MERTMEDICAL**  
**ENDOTEK®**

# Segue™ Multipurpose TransProtection System™

The Segue Multipurpose TransProtection System\* is designed for the transportation and protection of endoscopes and other critical devices pre- and post-procedure.



## Benefits:

- Complies with SGNA Standards for Infection Control\*\*
- Multipurpose design for use as a work surface, receptacle, and transportation system
- DriPAD™ lining traps fluid, keeping the surface dry
- Pad designed to absorb over 5 liters of fluid
- Transparency allows quick identification of contents
- Integrated peel-off label identifies Clean or Biohazard



\*Patent Pending

\*\*Standards of Infection Control in Reprocessing of Flexible Gastrointestinal Endoscopes, SGNA (c) 2012



## ORDERING INFORMATION

| Part Number | Description                         | Packaged |
|-------------|-------------------------------------|----------|
| SEG-1519    | Multipurpose TransProtection System | 5 Units  |