


INFLATION DEVICES

Simplifying Balloon Inflation



EXPANDING THE POSSIBILITIES



MERTMEDICA
ENDOTEK™

INFLATION DEVICES

COMPACT, EASY TO USE, ACCURATE

Merit Medical has achieved world market leadership in the inflation device market by listening and responding to customer needs. By pairing our inflation device experience with our fully covered non-vascular stent technology, Merit Endotek continues to find ways to meet the changing needs of our customers.



**ENDOTEK ANALOG
INFLATION DEVICE**
20ml , 30 ATM Inflation Device
with Angled Gauge

Catalog Number
ENDO-AN2030
Packaged 5 Per Box



**ENDOTEK DIGITAL
INFLATION DEVICE**
20ml , 30 ATM Inflation Device
with Backlit Gauge

Catalog Number
ENDO-DI2030
Packaged 5 Per Box



**ENDOTEK DIGITAL
ADVANCED INFLATION DEVICE**
20ml, 30 ATM Inflation Device
with New Enhanced Gauge Design

COMING SOON

INFLATION DEVICE FEATURES:

- **LOCK RELEASE HANDLE**
Allows for easy and accurate one handed control of inflation and deflation
- **HIGH PRESSURE FLEXIBLE TUBING**
Enables flexibility in positioning the device when attached to a balloon catheter
- **ANGLED PRESSURE GAUGE**
Designed to provide accurate pressure readings in low light environment
Available in both analog and digital configurations

OUR PRODUCT LINE



AEROSIZER™
Stent Sizing
Device



Fully Covered
ALIMAXX-ES™
Esophageal Stent



MAXXWIRE™
Guide Wire



Fully Covered
AERO™
Tracheobronchial
Stent



ENDOTEK Negative
Pressure Syringe



Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

www.merit.com

Merit Medical EUROPE, MIDDLE EAST AND AFRICA (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • The Netherlands • Tel: +31 43 358 82 40

Free phone for specific country:

Austria 0800 295374 • Belgium 0800 729 06 (Dutch) 0800 731 72 (French) • Denmark 80 88 00 24 • France 0800 91 60 30 • Germany 0800 182 08 71

Ireland (republic) 1800 553 163 • Luxembourg 8002 25 22 • Netherlands 0800 022 8184 • Sweden 020 79 24 45 • UK 0800 973 115 • Finland 00800 6374 8633

401887001/A ID 040710