

TWISTER[®] PLUS

Rotatable Retrieval Device

 **MERTMEDICAL**
ENDOTEK[®]

NEW: 15mm
PULMONARY
SPECIFIC SIZE

High Volume Net

360° Rotation

Kink-Resistant Catheter



TWISTER[®] PLUS

Rotatable Retrieval Device

- Pulmonary specific foreign body retrieval net
- Works through a bronchoscope 2.0mm working channel
- Larger sizes also available

Order Number	Description	Loop Diameter (mm)	Catheter OD (mm)	Working Length (cm)	Min. Working Channel (mm)	Packaging
TWIST15-PULM	TWISTER PLUS Rotatable Retrieval Device	15	1.7	120	2.0	5/Box

Merit Medical Endotek

1600 West Merit Parkway
South Jordan, Utah 84095
1-801-253-1600
1-800-35-MERIT
endotek.merit.com



EXPANDING THE POSSIBILITIES[™]



Please consult the product labels and inserts for any indications, contraindications, potential complications, warnings, precautions and directions for use.

TWISTER[®] PLUS is a trademark of Horizons International Corp. Protected by one or more of the following U.S. Patents: 9,101,342; 8,858,567

402603001/B ID 072117