

TWISTER[®] PLUS

Rotatable Retrieval Device

Available in
15mm, 22mm,
26mm and 32mm.

15mm size for use within
a 2.0mm working channel.

 **MERITMEDICAL**
ENDOTEK[®]

High Volume Net

360° Rotation

Kink-Resistant Catheter



TWISTER[®] PLUS

Rotatable Retrieval Device

NEW

Order Number	Description	Loop Diameter (mm)	Catheter OD (mm)	Working Length (cm)	Min. Working Channel (mm)	Packaging
TWIST15-PED	TWISTER PLUS Rotatable Retrieval Device	15	1.7	160	2.0	5/Box
TWIST22	TWISTER PLUS Rotatable Retrieval Device	22	2.5	230	2.8	5/Box
TWIST26	TWISTER PLUS Rotatable Retrieval Device	26	2.5	230	2.8	5/Box
TWIST32	TWISTER PLUS Rotatable Retrieval Device	32	2.5	230	2.8	5/Box

Merit Medical Endotek

1600 West Merit Parkway
South Jordan, Utah 84095
1-801-253-1600
1-800-35-MERIT
endotek.merit.com



EXPANDING THE POSSIBILITIES[™]



Please consult the product labels and inserts for any indications, contraindications, potential complications, warnings, precautions and directions for use.

TWISTER[®] PLUS is a trademark of Horizons International Corp. Protected by one or more of the following U.S. Patents: 9,101,342; 8,858,567

402604001/B ID 071417