

TIOTM THREE IN ONE

ORAL AIRWAY • BITE BLOCK • O₂ PORT



EXPANDING THE POSSIBILITIESTM

 **MERTMEDICAL**
ENDOTEK[®]

DO MORE WITH LESS.

The TIO (Three-In-One) is a single, elegantly-designed device that incorporates the benefits of an oropharyngeal airway (oral airway), bite block and oxygen administration port into one convenient device. The innovative design of the TIO enhances procedure efficiency and improves patient care during TEE and endoscopic procedures that require deep conscious sedation, and/or for patients with oropharyngeal obstructions such as sleep apnea.

WE KNOW WHAT YOU'RE THINKING: "WHY DIDN'T I THINK OF THIS?"






■ **UNIQUE ORAL AIRWAY DESIGN**
 offered in 90, 100 and 110mm blade lengths to meet anatomical requirements.




■ **COMFORTABLE 60F BITE BLOCK**
 designed to allow easy, safe insertion of a full range of endoscopes, flexible bronchoscopes, and other through-the-mouth instruments.




■ **CONVENIENT AND EFFECTIVE STANDARD OXYGEN PORT/CONNECTOR** designed to deliver supplemental oxygen up to 6.0L/min.*
*based on Merit Medical bench testing data

TIO THREE IN ONE		DISPOSABLE	
Catalog Number	Specifications (Packaged 10 units)		
	Bite Block Opening (F)	Oral Airway Length (mm)	Color
TIO60-90	60	90	
TIO60-100	60	100	
TIO60-110	60	110	


ALSO AVAILABLE



EndoMAXX® Esophageal Stent



AERO® Tracheobronchial Stent



ALIMAXX-B® Biliary Stent



BIG60® INFLATION DEVICE



endotek.merit.com

