


ELATION™

BALLOON DILATOR

Single-Use Fixed Wire Balloon Dilator

Merit Medical Endotek's Elation Fixed Wire Balloon Dilators provide precise, controlled, three-stage dilation and crystal clear visualization of tight strictures during balloon endoscopy.



● **Unique Catheter Design**

Provides excellent pushability, trackability and rapid deflation

● **Optimal Balloon Clarity**

For exceptional endoscopic visibility of the stricture during dilation

● **Radiopaque Marker Bands**

Allow fluoroscopic identification during balloon placement

● **Proprietary Balloon Design**

Ensures three specific diameters at three corresponding pressures, providing accurate, reliable and repeatable balloon sizing



Single-Use Fixed Wire Balloon Dilator

Specifications (Packaged 1 Unit)

Catalog Number	Inflated O.D. (mm)	Inflated O.D. (F)	Balloon Length (cm)	Inflation Pressures (ATM)	Inflation Pressures (kPa)	Min. Working Channel (mm)	Catheter Size (F)	Working Length (cm)
E6	6-7-8	18-21-24	8	6-8.5-10	608-861-1013	2.8	7.5F	180
E8	8-9-10	24-27-30	8	4-6-8	405-608-811	2.8	7.5F	180
E10	10-11-12	30-33-36	8	4.5-6.5-8	456-659-811	2.8	7.5F	180
E12	12-13.5-15	36-40.5-45	8	4-6.5-8	405-659-811	2.8	7.5F	180
E15	15-16.5-18	45-49.5-54	8	4.5-5.7	405-557-709	2.8	7.5F	180
E18	18-19-20	54-57-60	8	3.5-5-6	355-507-608	2.8	7.5F	180



EXPANDING THE POSSIBILITIES™

Merit Medical Endotek

1600 West Merit Parkway
South Jordan, Utah 84095
1-801-253-1600
1-800-35-MERIT
endotek.merit.com

Merit Medical Europe, Middle East, and Africa (EMEA)

Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
Tel: +31 43 358 82 22

Merit Medical Ireland Ltd.

Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Free phone for specific country

Austria	0800 295 374	Luxembourg	8002 25 22
Belgium	0800 72 906 (Dutch)	Netherlands	8000 022 81 84
	0800 73 172 (French)	Norway	800 11629
Denmark	80 88 00 24	Sweden	020 792 445
France	0800 91 60 30	Switzerland	800 9000 5060
Finland	0800 770 586	UK	0800 973 115
Germany	0800 182 0871		
Ireland (Republic)	1800 553 163		
Italy	800 897 005		



403153002/A ID 011215