


# ENDOTEK **BIG60**<sup>TM</sup>

Simplifying  
Balloon  
Dilation



EXPANDING THE POSSIBILITIES<sup>TM</sup>



**MERTMEDICA**  
**ENDOTEK<sup>®</sup>**

# ENDOTEK **BIG60**<sup>TM</sup> INFLATION DEVICE

*Merit Endotek's BIG60 Inflation Device is designed to inflate and deflate non-vascular balloon dilators while accurately monitoring and displaying crucial inflation pressures. Our ergonomically designed 60ml inflation device is a one-piece, easy-to-use device that pairs Merit's market-leading inflation device experience with our non-vascular stent technology.*



**ENDOTEK BIG60 INFLATION DEVICE**  
60ml Inflation Device with Angled Gauge

Catalog Number: **ENDO-AN6012**  
Packaged 5 Per Box



## INFLATION DEVICE FEATURES:

- **CONVENIENT ALL-IN-ONE DESIGN**  
*Requires no assembly or cleaning before use and is completely disposable*
- **LOCK RELEASE HANDLE**  
*Allows for easy and accurate one-handed control of inflation and deflation*
- **FLEXIBLE HIGH PRESSURE TUBING**  
*Aids in priming and enables flexibility in positioning the device when attached to a balloon catheter*
- **ANGLED ANALOG PRESSURE GAUGE**  
*Provides accurate pressure readings and high visibility in low-light environments*

### ALSO AVAILABLE



Endotek Analog Inflation Device  
20ml, 30ATM



Endotek Digital Inflation Device  
20ml, 30ATM



Brighton<sup>®</sup> Bipolar Coagulation Probe



Fully Covered ALIMAXX-ES<sup>™</sup> Esophageal Stent



ALIMAXX-B<sup>®</sup> Biliary Stent



Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

www.merit.com

Merit Medical EUROPE, MIDDLE EAST AND AFRICA (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • The Netherlands • Tel: +31 43 358 82 40

Free phone for specific country:

Austria 0800 295374 • Belgium 0800 729 06 (Dutch) 0800 731 72 (French) • Denmark 80 88 00 24 • France 0800 91 60 30 • Germany 0800 182 08 71

Republic of Ireland 1800 553 163 • Luxembourg 8002 25 22 • Netherlands 0800 022 8184 • Sweden 020 79 24 45 • UK 0800 973 115 • Finland 00800 6374 8633

402007001/B ID 050411