


VACLOK®

Negative Pressure Syringe



EXPANDING THE POSSIBILITIES



MERITMEDICA
ENDOTEK™


VACLOK®

Negative Pressure Syringe


Merit Endotek's VacLok Negative Pressure Syringe is the perfect tool for creating negative vacuum pressure. With the VacLok, Merit Endotek continues to support clinicians by turning negative pressure into a positive experience.



Step 1: Pull



Step 2: Rotate & Release



ALIMAXX-ES®
Fully Covered
Esophageal Stent



MAXXWIRE®
Guide Wire



Brighton® Bipolar
Coagulation Probe



BIG60™ Inflation Device



AERO® Fully Covered
Tracheobronchial Stent

SYRINGE FEATURES:

HOLDS NEGATIVE PRESSURE

Plunger design holds negative pressure with ease, and won't inadvertently disengage

ERGONOMIC DESIGN

Large wings and skid resistant thumb pad make syringe easy to grip and control

VARIOUS SIZES AVAILABLE

Available in 20ml, 30ml, 60ml

DUAL USE DEVICE

Can be used as negative pressure syringe or as standard syringe

NEGATIVE PRESSURE SETTINGS

Four or six different negative pressure settings (dependent on syringe size)

VacLok Negative Pressure Syringe

Catalog Number	Size	Pres. Settings	Package
ENDO-VAC20	20ml	4 position	20/box
ENDO-VAC30	30ml	6 position	20/box
ENDO-VAC60	60ml	6 position	20/box



Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

www.merit.com/endotek

Merit Medical EUROPE, AFRICA AND THE MIDDLE EAST (EMEA) Amerikalaan 42, 6199 AE Maastricht-Airport, The Netherlands Tel: +31 43 358 82 22

Merit Medical Ireland Ltd. • Parkmore Business Park West Galway, Ireland • +353 (0) 91 703 733

Free phone for specific country: Austria 0800 295 374 • Belgium 0800 72 906 (Dutch) 0800 73 172 (French) • Denmark 80 88 00 24 • France 0800 91 60 30

Germany 0800 182 0871 • Ireland (republic) 1800 553 163 • Luxembourg 8002 25 22 • Netherlands 0800 022 81 84 • Sweden 020 792 445 • UK 0800 973 115

402284001/A ID 090711